




E. Coli as a host organism in molecular biology




Electron micrograph of an E. coli cell ruptured to release its DNA. The tangle is a portion of a single DNA molecule containing over 4.6 million base pairs encoding approximately 4,300 genes. The small circlets are plasmids.

(Courtesy of Huntington Potter and David Dressler, Harvard Medical School, accessed at <http://users.rcn.com/jkimball.ma.ultranet/BiologyPages/R/RecombinantDNA.html>)




Plasmid DNA: extrachromosomal circular molecules in the bacterial cell



<http://w3.uniroma1.it/centric-nr-can/edocs/gme3.htm>

relaxed

supercoiled



<http://intranet.siu.edu/~mbmb/451resources/supercoiling.html>

0.2 μm


Shigella Outbreaks in Japan, 1953-1960: occurrence of antibiotic-resistant strains

Year	# Strains tested	Sm	Tc	Cm	Sm + Cm	Sm + Tc	Cm + Tc	Sm, Cm & Tc
1953	4900	5	2	0	0	0	0	0
1954	4876	11	0	0	0	0	0	0
1955	5327	4	0	0	0	0	0	1
1956	4399	8	4	0	0	0	1	0
1957	4873	13	46	0	2	2	0	37
1958	6563	18	20	0	7	2	0	193
1959	4071	16	32	0	71	0	0	74
1960	3396	29	36	0	61	9	7	308

(From <http://www.mun.ca/biochem/courses/4103/topics/plasmids.html>)


- Plasmids confer antibiotic resistance to strains harboring them!
- Plasmids can be “shared” between strains by lateral transfer!

Another View: Plasmid Maps



This plasmid is a bacterial:yeast shuttle vector!

Plasmids can be used as vectors to carry a gene of interest




Cloning into a plasmid

(from <http://www.accessexcellence.org/RC/VL/GG/plasmid.html>)

Gene Cloning I

- Provides an easily replicated source of a gene.
- Enables detailed analysis of gene structure.
- Permits production of large quantities of the gene product for functional and structural analysis.




Courtesy of U.S. Department of Energy Genome Program




Taking a closer look at our plasmid, pMSH2

MSH2 construct for mutagenesis and functional assays




-produces a protein product of the expected molecular weight




-complements *msh2* Δ defects, therefore the construct is functional

Strategy for subcloning DNA fragments into a plasmid vector



Inserting a DNA Sample into a Plasmid


Another View: Plasmid Maps



What are restriction sites for??

Biological function of restriction endonucleases is to protect cells from foreign DNA

Infecting DNA is cleaved (restricted) by the restriction enzymes, preventing it from successfully replicating and parasitizing the cell



most lab strains are completely “domesticated” (R-M systems have been inactivated)

some examples in the genotypes:


hsd

mcrA

mcrBC

Restriction Endonucleases


- Over 10,000 bacteria and archaea have been screened for restriction enzymes
- Restriction enzymes are not confined exclusively to bacteria
- Nearly 3,000 enzymes have been found, exhibiting over 200 different specificities (many of the 3000 are isoschizomers; different enzyme, same recognition site)



***EcoRI* binding DNA**

http://www.biophysics.pitt.edu/john_m_rosenberg.htm

Activity of a typical Restriction Enzyme




**Restriction Enzyme
Action of EcoRI**

ends generated from a Type II restriction endonucleases

	enzyme	recognition site	end generated		
defined ends	<i>PvuII</i>	5' CAGCTG 3' GTCGAC	5' CAG 3' GTC	CTG GAC	blunt
	<i>BamHI</i>	5' GGATCC 3' CCTAGG	5' G 3' CCTAG	GATCC G	5' overhang
	<i>PstI</i>	5' CTGCAG 3' GACGTC	5' CTGCA 3' G	G ACGTC	3' overhang
variable ends	<i>XmnI</i>	5' GAANNNTTC 3' CTTNNNAAG	5' GAANN 3' CTTNN	NNTTC NNAAG	blunt
	<i>BanI</i>	5' GGPuPuCC 3' CCPuPyGG	5' G 3' CCPuPyG	GPuPuCC G	5' overhang
	<i>BstXI</i>	5' CCANNNNNTGG 3' GGTNNNNNACC	5' CCANNNNN 3' GGTN	NTGG NNNNNACC	3' overhang

Restriction Enzymes I

- Restriction enzymes can be used:
 - to cut DNA at specific nucleotide sequences.
 - Example – to cut chromosomes into smaller pieces for analysis by gel electrophoresis.
 - for cloning or the generation of genetic libraries.




Courtesy of U. S. Department of Energy



Restriction Enzymes II


- Can be used to detect mutations in DNA.
 - Example – The enzyme called DdeI can identify the mutation that causes sickle cell anemia.
 - The mutation changes the DNA sequence so that DdeI cannot cut the DNA if the mutation is present.




Courtesy of Alford, Rossiter, and Caskey.



We will use restriction digestion to distinguish between WT and mutagenized alleles of our pMSH2 vector




After the diagnostic gel, only proceed with a plasmid showing the altered restriction endonuclease digestion pattern (e.g. 2 or 5)



✓ Inoculate 25 ml of media with the correct bacterial colony for a midi-scale plasmid preparation (anion exchange chromatography)

✓ Use miniprep DNA to transform yeast to begin the functional analyses

Project Overview



Manipulate the yeast *MSH2* gene to determine which human missense mutations are likely to be benign or pathogenic in nature.

Examine the defect at a molecular level to determine why the Msh2 variants are dysfunctional.

Locating the codon to be mutagenized

Note that the two sequences don't have corresponding codon numbers

for example Ala 719 in yeast is Ala 700 in humans

Yeast	600	INITLTYTPVFEKLSLVLAHLDVIA SFAHTSSYAPIPYIRPKLHPMDSERRTHLISSRHP	659
		+NI+ Y + L+ VLA LD + SFAH S+ AP+PY+RP + + R L +SRH	
Human	582	VNISSGYVEPMQTLNDVLAQLDAVVSFAHVSNGAPVPYVRPAILE-KGQGRIILKASRHA	640
Yeast	660	VLEMQDDISFISNDVTLESGKGDFLIITGPNMGGKSTYIRQVGVISLMAQIGCFVPCEE	A 719
		+E+QD+I+FI NDV E K F IITGPNMGGKSTYIRQ GVI LMAQIGCFVPCE	A
Human	641	CVEVQDEIAFIPNDVYFEKDKQMFHIIITGPNMGGKSTYIRQTGVIVLMAQIGCFVPCES	A 700
Yeast	720	EIAIVDAILCRVGAGDSQLKGVSTFMVEILETASILKNASKNSLIIVDELGRGTSTYDGF	779
		E++IVD IL RVGAGDSQLKGVSTFM E+LETASIL++A+K+SLII+DELGRGTSTYDGF	
Human	701	EVSIVDCILARVGAGDSQLKGVSTFMAEMLETASILRSATKDSLIIIDELGRGTSTYDGF	760
Yeast	780	GLAWAIAEHIAASKIGCFALFATHFHELTELSEKLPNVKNMHVVAHIEKNLKEQKHDDEDI	839
		GLAWAI+E+IA+KIG F +FATHFHELT L+ ++P V N+HV A +E +	
Human	761	GLAWAISEYIATKIGAFCMFATHFHELTALANQIPTVNNLHVTA-----LTTEETL	811
Yeast	840	TLLYKVEPGISDQSFGIHVAEVVQFPEKIVKMAKRKANELDDLKTNNE----DLKK--AK	893
		T+LY+V+ G+ DQSFGIHVAE+ FP+ +++ AK+KA EL++ + E D+ + AK	
Human	812	TMLYQVKKGVCDQSFGIHVAELANFPKHVIECAKQKALELEEFQYIGESQGYDIMEPAAK	871

Designing a mutagenic oligonucleotide - choosing the codon guided by a yeast codon usage table

yeast codon 122 changed from **glycine (G)** to **serine (S)**
termed G122S mutation

wild-type *MSH2* coding sequence

5' ^ACGCATCTCCA^GGAACATTGAGC
^S ^P ^N ^I ^E



missense mutation in *MSH2* coding sequence

5' ^ACGCATCTCCA^STCTAACATTGAGC
^S ^P ^N ^I ^E

5' ^ACGCATCTCCA^STCAAACATTGAGC
^S ^P ^N ^I ^E

5' ^ACGCATCTCCA^STCCAACATTGAGC
^S ^P ^N ^I ^E

5' ^ACGCATCTCCA^SAGTAACATTGAGC
^S ^P ^N ^I ^E

top four candidates

Codon usage in yeast (Appendix of manual)

AmAcid	Codon	Number	/1000	Fraction
Ser	AGT	2411.00	11.14	0.14
Ser	AGC	1599.00	7.39	0.09
Ser	TCG	1411.00	6.52	0.08
Ser	TCA	3316.00	15.33	0.19
Ser	TCT	5495.00	25.40	0.32
Ser	TCC	3198.00	14.78	0.18

Designing a mutagenic oligonucleotide:

choosing the codon to create or destroy a restriction endonuclease site

wild-type *MSH2* coding sequence

```
      A   S   P   G   N   I   E  
5' CGCATCTCCAGGGAACATTGAGC 5'  
   GCGTAGAGGTCCCTTGTAACTCG
```



lose a *Bst*NI site

missense mutation in *MSH2* coding sequence

```
      A   S   P   S   N   I   E  
5' CGCATCTCCATCTAACATTGAGC 5'  
   GCGTAGAGGTAGATTGTAACTCG
```

✓ best codon choice

✓ lose a *Bst*NI site

Designing a mutagenic oligonucleotide - choosing the codon guided by a yeast codon usage table

yeast codon 122 changed from **glycine (G)** to **serine (S)**
termed G122S mutation

wild-type *MSH2* coding sequence

5' ^ACGCATCTCCA^GGAACATTGAGC^N^I^E



missense mutation in *MSH2* coding sequence

5' ^ACGCATCTCCA^STCTAACATTGAGC^N^I^E

5' ^ACGCATCTCCA^STCAAACATTGAGC^N^I^E

5' ^ACGCATCTCCA^STCCAACATTGAGC^N^I^E

5' ^ACGCATCTCCA^SAGTAACATTGAGC^N^I^E

top four candidates

Codon usage in yeast (Appendix of manual)


AmAcid	Codon	Number	/1000	Fraction
Ser	AGT	2411.00	11.14	0.14
Ser	AGC	1599.00	7.39	0.09
Ser	TCG	1411.00	6.52	0.08
Ser	TCA	3316.00	15.33	0.19
Ser	TCT	5495.00	25.40	0.32
Ser	TCC	3198.00	14.78	0.18

Laboratory Session 3 (Mon/Tue)

designing the mutagenic oligonucleotide

important considerations:

- 23 nucleotide long (for efficient annealing to the template)



- codon preference for expression in yeast
- creating/destroying a restriction site for screening
- strand must be complementary to template