

PCR reaction

Chapter 20


Polymerase chain reaction (PCR)

- Copy DNA without vectors and host cells
- Goal to make many copies of DNA in a defined region
- Uses high concentration of two primers that are complementary to sequences at the ends of the DNA region to be amplified, deoxynucleoside triphosphates (dNTPs), and a heat-stable form of DNA polymerase called *Taq* polymerase
- Sample of DNA taken through repeated cycles of denaturation, annealing and synthesis
 - Thermocycler automates this process
- After 30 cycles of amplification, a DNA sample will increase 2³⁰-fold






Copyright © The McGraw-Hill Companies, Inc. Permission required for reproduction or display. 3' 5' Start with doublestranded template **Primer** DNA, present in annealing Forward primer low amounts, plus site 2 primers present in high amounts, dNTPs, and Taq polymerase. Cycle 1 **Denaturation: Heat the** DNA to separate strands. Primer annealing: Lower the temperature to allow primers to bind to the template DNA. Reverse primer **Primer extension: Incubate** at a higher temperature to allow the synthesis of the complementary strand.




With each successive cycle, the relative amount of DNA fragments that end exactly at both primer sites (marked *) increases. Therefore, after many cycles, the vast majority of DNA fragments contain only the region that is flanked by the 2 primer sites.

DNA sequencing

- Determines base sequence of DNA
- Dideoxy chain-termination method or dideoxy sequencing
 - Dideoxynucleoside triphosphates (ddNTPs) are missing the 3' –
 OH group and will terminate the chain
 - 4 tubes with many copies of single stranded DNA of interest
 - Each tube has a different radiolabelled dNTP
 - DNA polymerase will make complementary strand until dNTP inserted and chain terminates
 - After electrophoresis, DNA sequence can be read by reading which base is at the end of the DNA strand
- Procedure has been automated using fluorescent dyes in one tube






GENOMES & PROTEOMES

A Microarray Can Identify Which Genes Are Transcribed by a Cell


- DNA microarray or gene chip
- Used to monitor the expression of thousands of genes simultaneously
- Short sequences of known genes attached to spots on slide
- Goal to find out which genes are transcribed into mRNA in particular sample of cells
- mRNA isolated from those cells and used to make fluorescently labeled cDNA
- cDNAs that are complementary to the DNAs in the microarray will hybridize
- If the fluorescence intensity in a spot is high, a large amount of cDNA was in the sample that hybridized to the DNA at this location

GENOMES & PROTEOMES



DNA fingerprinting

- Identifies and distinguishes among individuals based on variations in their DNA
- Chromosomal DNA produces series of bands on a gel
- Unique pattern of bands used
- Automated using PCR to amplify short tandem repeat sequences (STRs) aka microsatelite DNA
 - Such tandem repeat sequences are found at specific locations in the genomes of all species, and the number of repeats at each spot tends to vary from one individual to the next



Uses

- Identify different species of bacteria and fungi
- Forensics 1986 first use in US court system
- Paternity testing and other family relationships